

826 – Blue:

The advantage of the 826 alginate is the long setting time of 8 minutes with a water temperature of 26° C, the operator is not rushed even for big projects like torso castings and the temperature (26° C) makes the model more comfortable than using cold water.

It sets firm and strong and has great details in the casting

Mixing ratio: 3,75-4 : 1 (water : powder, by weight)

Working time: **8 minutes**



Tips:

*) If you do not use your alginate casting immediately, put some wet towels into the negative because the alginate dries up very quickly, shrinks and is not useable anymore.

Shrinkage – (with the use of wet towels)

10 minutes	< 0,10%
20 minutes	approx. 0,16 %
30 minutes	approx. 0,33 %
40 minutes	approx. 0,46 %
50 minutes	approx. 0,62 %
60 minutes	approx. 0,75 %

*) As always - warmer water will accelerate the setting time, cooler water will retard the time.

Blond & Braun

Dr. Hans Lechner Str. 7, A – 5071 Wals

Tel.: +43 – (0)662 – 85 02 02 – 0, Fax: +43 – (0)662 – 85 02 02 – 11

mailto: office@blond-braun.at / web: www.blond-braun.at